

i.Mune NEO – A New Test For Epigenetic Immune Cell Quantification in Newborns

Schulze J.¹, Lundbäck, V.², Werner J.¹, Walter S.¹, Staub U.¹, Sachsenmaier C.¹

¹ Epimune GmbH, Berlin Germany


² Karolinska University Hospital, Centre for Inherited Metabolic Diseases

Introduction

Inborn Errors of Immunity (IEI) are characterized by quantitative dysregulation of immune cells already at birth. More than 400 IEIs have been described that affect almost all components of the immune system.

Current newborn screening for IEI is limited to quantification of T-cells using the **T-Cell Receptor Excision Circle (TREC)** test and to B-cells using the **Kappa Receptor Excision Circle (KREC)** test. Both tests produce a fair number of non-actionable screening findings and false-positive results.

Also, expanded screening for quantitative abnormalities of more immune cells has the potential to significantly improve early detection of more IEIs in newborns than currently possible.




©Roshini S. Abraham, PhD

Approach

We developed real-time PCR based **epigenetic immune cell quantification** from a drop of dried blood.

This offers **confirmatory testing** possibilities to existing screening methods as well as the potential to **expand newborn screening** to a variety of other IEIs.




Results

Successful identification of SCID and XLA cases in newborn DBS samples

Identifier	Classification	Gene Defect	Loss of Function Type	TREC/KREC Newborn Screening			Epigenetic qPCR Analysis			
				TREC ⁽¹⁾ Positive [yes/no]	KREC ⁽²⁾ Positive [yes/no]	Screening Classification	(CD3 G/D, GAPDH) ⁽³⁾ Conspicuous [yes/no]	(MVD, GAPDH) ⁽³⁾ Conspicuous [yes/no]	(LRPS, GAPDH) ⁽³⁾ Conspicuous [yes/no]	Screening Classification
1	SCID	ADA	amorph	yes	yes	correctly identified	yes	yes	yes	correctly identified
2	SCID	ADA	amorph	no	yes	correctly identified	yes	yes	yes	correctly identified
3	DO-SCID ⁽⁴⁾	ADA	hypomorph	no	yes	correctly identified	no	yes	yes	correctly identified
4	DO-SCID ⁽⁴⁾	ADA	hypomorph	no	yes	correctly identified	yes	yes	yes	correctly identified
5	SCID	AK2	amorph	yes	no	correctly identified	yes	yes	yes	correctly identified
6	SCID	AK2	amorph	yes	yes	correctly identified	yes	yes	no	correctly identified
7	SCID	Artemis	amorph	yes	yes	correctly identified	yes	yes	yes	correctly identified
8	SCID	CD3D	amorph	yes	no	correctly identified	yes	yes	no	correctly identified
9	SCID w ME ⁽⁵⁾	IL2RG	amorph	yes	no	correctly identified	no	no	no	not identified
10	SCID	IL2RG	amorph	yes	no	correctly identified	yes	yes	yes	correctly identified
11	SCID	IL7RA	amorph	yes	no	correctly identified	yes	no	no	correctly identified
12	SCID	IL7RA	amorph	yes	no	correctly identified	yes	yes	yes	correctly identified
13	SCID	IL7RA	amorph	yes	no	correctly identified	yes	yes	yes	correctly identified
14	DO-SCID ⁽⁴⁾	JAK3	hypomorph	no	no	not identified	yes	yes	yes	correctly identified
15	SCID	PNP	amorph	yes	yes	correctly identified	yes	yes	yes	correctly identified
16	SCID	PNP	amorph	yes	yes	correctly identified	yes	yes	yes	correctly identified
17	SCID	RAG1	hypomorph	yes	yes	correctly identified	yes	yes	no	correctly identified
18	SCID	RAG1	yes	yes	yes	correctly identified	no	yes	yes	correctly identified
19	SCID	RAG2	amorph	yes	yes	correctly identified	yes	no	yes	correctly identified
20	XLA	BTX	amorph	no	yes	correctly identified	yes	no	yes	correctly identified
21	XLA	BTX	amorph	no	yes	correctly identified	no	no	yes	correctly identified
22	XLA	BTX	amorph	no	yes	correctly identified	no	no	yes	correctly identified
23	XLA	BTX	amorph	no	yes	correctly identified	yes	yes	yes	correctly identified
24	XLA	BTX	hypomorph	no	no	not identified	no	yes	yes	correctly identified

Baron et al., 2018

Reduction of false-positive TREC screening referrals



Blom et al., 2021


i.Mune NEO - a real-time PCR based test for immune cell quantification in dried blood spot samples from newborns

	i.Mune NEO	i.Mune NEO LITE	i.Mune NEO CD3
T lymphocytes (CD3+)	✓	✓	✓
Memory Helper T lymphocytes (CD3+CD4+CD45RO+)	✓	✓	✓
B lymphocytes (CD19+)	✓	✓	
Natural Killer lymphocytes (CD16+CD56dim)	✓	✓	
Helper T lymphocytes (CD3+CD4+)	✓		
Cytotoxic T lymphocytes (CD3+CD8+)	✓		

Summary & Outlook

i.Mune NEO is a novel epigenetic immune cell quantification test allowing confirmatory testing of TREC- and KREC positive screening cases. It also is suitable for screening applications (i.Mune NEO LITE, NEO CD3). Going forward, addition of other epigenetic immune cell markers will allow expanded screening for additional IEI – e.g. Severe Congenital Neutropenia (SCN), Immune Dysregulatory Disorders, NK deficiencies and others (see also: Blom et al., Int. J. Neonatal Screen. 2021, 7, 74)

Ongoing evaluation of i.Mune NEO in a cohort of newborn screening samples (Karolinska University Hospital)



Lundbäck (unpublished)